

Product Datasheet

www.leinco.com



**Leinco
Technologies, Inc.**

excellence in early discovery research

8X ClearFlow™ Sheath Fluid

“a cleaner stream”

Prod. No.: S622
Lot No.: **Always Refer to Lot Specific Datasheet**
Pkg. Size: 2 x 2.5 gallons

Description

Background:

8X ClearFlow™ Sheath Fluid is specifically formulated for transporting particles through a flow cytometer. ClearFlow™ is specially manufactured in our class 10,000 clean room environment for low particule count and fluorescence backgrounds to guarantee better signal to noise ratio measurements. **This reagent is for *in vitro* diagnostic use only.**

Recommended Use:

This product contains an antibacterial/antifungal preservative (other than sodium azide) and may not be suitable for cell sorters and systems requiring a preservative free buffer.

Direction for Reconstitution and Use:

PLEASE REFER TO LOT SPECIFIC DATA SHEET FOR FILLING INTSTUCTIONS

Note: To avoid issues using new sheath brands sheath filters need to be cleaned or changed and then system flushed with DI Water for approximately 30 minutes.

Storage and Stability

Store at room temperature. **Do not Freeze.**

Country of Origin

USA



Products are for research use only. Not for use in diagnostic or therapeutic procedures.