## **Product Data Sheet Native Antigen**



| Product name          | Chlamydia trachomatis Antigen   |  |
|-----------------------|---|--|
| Order no.             | BA1372VSER  |  |
| Lot no.               | Lot specific  |  |
| Agent                 | Chlamydia trachomatis, Lymphogranuloma venereum (LGV II)<br>strain 434  |  |
| Culture system        | Vero cell culture   |  |
| Preparation           | Chlamydia trachomatis elementary and reticular bodies, purified from infected host cells by density gradient centrifugation.  Presented in phosphate buffered saline, pH 7.5.   |  |
| Inactivation          | Heat treatment and sonication. Inactivity is confirmed by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions. |  |
| Preservatives         | None  |  |
| Protein concentration | Lot specific,<br>specified in mg/mL   | Determined by Bradford assay (BioRad) using BSA as standard. |
| Analysis              | By SDS-PAGE and CFT   |  |
| Storage               | -65 °C or lower   | Avoid repeated freeze-thaw cycles. Sonicate before use.      |
| Date of manufacture   | Lot specific  |  |
| Expiry date           | 10 years from date of manufacture (under recommended storage conditions)  |  |
| Manufacturer          | Institut Virion\Serion GmbH, Germany  |  |
|                       |   |  |

- ✓ For *in vitro* investigational use only. Not for use in diagnostic or therapeutic procedures.
- ✓ If skin contact occurs, wash contacted area with soap and water. If irritation occurs, consult a physician.