

The right solution for guideline compliant practice

SERION ELISA *control*

SERION ELISA *controls* are positive control sera for the qualitative and quantitative antibody determination when using immunoassays of the SERION ELISA *classic* and SERION ELISA *antigen* product line. The ready-to-use components are, in addition to the controls supplied with the SERION ELISA test kits, supplementary external reagents to determine validity of SERION ELISA test runs, precision and reliability of the method.

In comprehensive validation studies, target values and validity ranges for each test-specific SERION ELISA *control* are determined and documented on lot-specific certificates of analysis. These values are recommended as references for the quality management in accredited laboratories when using SERION ELISA *controls* in combination with SERION ELISA test kits.

SERION ELISA *controls* are ready-to-use and must not be diluted any further. Rf-absorption must also not be performed.

Highlights

- Prediluted reagents allow easy handling
- Additional quality assurance through determination of methodical precision
- Laboratory internal error limits can be determined for an optimized quality control performance



Downloads

Flyer

Instructions for Use

QC Certificate

Safety Data Sheet

serion\diagnostics

[Imprint](#)
[Privacy policy](#)

[Virion\Serion Website](#)
[SERION Immunologics Website](#)

Contact us

[National Sales](#)
[International Sales](#)
[Customer Support](#)
[Product management](#)

Distributor Access

[Registration](#)
[Login](#)