

# **Neutrophil Elastase Monoclonal Antibodies**

# **ORDERING INFORMATION**

Catalog No.	Clone No.	MAb Subtype	Size	Library Pack No.	100ug/clone
13201	ELA11-207.2	IgM <sup>1</sup>	100ug, 500ug	132100	All 4 clones
13202	ELA11-207.3	IgM <sup>1</sup>	100ug, 500ug		
13203	ELA10-101.5	lgG1 <sup>2</sup>	100ug, 500ug		
13204	ELA10-103.2	lgG2b <sup>2</sup>	100ug, 500ug		
<sup>1</sup> PEG-purified antibody in PBS, pH 7.4					
<sup>2</sup> Protein G-purified antibody in PBS, pH 7.4					

# BACKGROUND

Neutrophil elastase is a major neutral protease present in human neutrophils and is released by neutrophils following stimulation with a variety of stimuli. It is believed to be a major mediator of the lung destruction seen in emphysema. In addition to elastin, neutrophil elastase can cleave other biologically active proteins such as fibrinogen and antithrombin and is released by neutrophils during clot formation.

# SPECIFICATION SUMMARY

Antigen: Purified human neutrophil elastase.

### Host Species: Mouse

**Specificity:** These antibodies recognize human neutrophil elastase and vary in their ability to inhibit neutrophil elastase activity as determined in an assay of elastase-mediated hydrolysis of elastin-Congo Red substrate:

Catalog No.	% Elastase Inhibition		
13201	1.2		
13202	4.1		
13203	90.0		
13204	82.0		

# APPLICATIONS

These antibodies have been qualified for use in ELISA to detect human neutrophil elastase and in studies of inhibition and enzyme kinetics of human neutrophil elastase.

### **DILUTION INSTRUCTIONS**

Dilute in PBS or medium that is identical to that used in the assay system.

# STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.