

Bovine Osteocalcin Monoclonal Antibody

ORDERING INFORMATION

Catalog No.: 50120 Clone No.: OCG4 Size: 100ug Format: Antibody was purified by column chromatography, dissolved in 10mM PBS, pH7.4, containing 1.0% bovine serum albumin, and then lyophilized. The lyophilized antibody does not contain preservative.

BACKGROUND

Osteocalcin is a noncollagenous, 49 amino acid glutamate-rich polypeptide bone matrix protein with a molecular weight of about 5800Da. Osteoblasts produce osteocalcin and incorporate it into the bone matrix. Osteocalcin is released into the circulation from the matrix during bone resorption and, therefore, is considered a marker of bone turnover rather than a specific marker of bone formation.

SPECIFICATION SUMMARY

Antigen: Bovine osteocalcin

Host Species: Mouse

Antibody Subtype: IgG1

Specificity: Epitope within residues 4-9; human, bovine, rabbit, pig, dog, goat, sheep, and chicken.

APPLICATIONS

Sandwich ELISA: 1-10ug/ml with immobilized antigen, pair with catalog no. 50116. Immunoblotting: 1-10ug/ml, reducing or non-reducing conditions.

Immunohistochemistry: 1-10ug/ml, frozen or paraffin-embedded tissue after proteinase K (0.4m/ml) antigen retrieval.

Chondroma Synovial Sarcoma Synovial Sarcoma




Figure. Immunohistochemistry: 1-10ug/ml, frozen or paraffin-embedded tissue after proteinase K (0.4m/ml) antigen retrieval. Chondroma: Human bone tumor. Synovial Sarcoma: human osteosarcoma.

Endusers should determine optimal concentrations for their applications.

DILUTION INSTRUCTIONS

Dissolve lyophilized antibody in 50ul of distilled water (2.0mg/ml final concentration) then dilute further for applications in PBS or media that is identical to that used in the assay system.

STORAGE AND STABILITY

Store stock solution (2.0 mg/ml) in aliquots at -20°C for 1 year, or at 4°C for 6 months after adding 0.1% sodium azide. Avoid repeated freeze-thaw cycles. Diluted antibody should not be stored.

Products are for research use only. Not for use in diagnostic or therapeutic procedures. QED Bioscience, Inc. 410 Axminister Dr, St. Louis, MO 63026 Phone: (800) 538-1145 Email: leincoglobal@leinco.com